

COVIDUR FACE SHIELD USER'S MANUAL



- **Covidur Face Shield** is a standalone device used to protect the user's eyes and face from bodily fluids, liquid splashes, or potentially infectious materials. Covidur face shield is situated at the crown of the head and is constructed with plastic to cover the user's **Eyes , Nose , Mouth and Face**.
- **Covidur Face Shield** can be intended for Medical purposes. Covidur Face Shield is effective at preventing HCP exposure to fluid biological airborne particulates.
- **Covidur Face Shield** is a single user equipment and can be single use, or can be multiple uses by the same user.
- **Covidur Face Shield** can be used easily in Daily Life, and it can be washed and disinfected with suitable product . Instructions for recommended cleaning and/or disinfection materials and processes are below
- It can be hand washed, disinfected using a 50% alcohol disinfectant.
- **Covidur Face Shield** does not contain any materials that will cause flammability, (not recommended against use in the presence of high intensity heat source or flammable gas)
- **Covidur Face Shield** does not alone will prevent infection from microbes or viruses. As indicated in Section 2, face shields may be effective at preventing HCP exposure to certain particulates during face shield shortages by providing barrier HCP protection to the wearer during COVID-19.
- While all the other visors are open with the top and sides of head, **Covidur Face Shield** is closed except chin. This provides complete protection.
- The feature of **Covidur Face Shiled** is that glasses and visor are together and it has an area that can be breathable easily from the mouth and nose.
- **Covidur Face Shiled** is unbreakable and shatter proof structure, it can be used for a long time.
- **Covidur Face Shiled** is certified from international organizations(**ISO 13485:2016** and **CE 93/42/EEC MEDICAL DEVICE DRECTIVE**). It is produced in an institution authorized by Republic of Turkey Ministry of Health.
- **Covidur Face Shiled** contains anti-fog material. It dries itself when breathing through the nose.


MATERIALS USED IN THE PRODUCT

Polystyrene.

ATTENTION OR CAUTIOUS WHILE USING

- It has sharp edges and corners, it can be easily scratched. It is recommended to be kept clean and dry and not to be exposed to direct sunlight when not in use.
- It is not recommended to use while driving.
- It can cause headaches and discomfort with your eyes during long-term use. It is suitable for short-term wear&take off usage.
- It should be used with caution for those with ocular problems and those who wear glasses. Those who use glasses can wear their glasses from outside the mask.
- It is not recommended for children.
- Keep away from chemicals and scratchers.
- Evaporation may occur during long-term use.
- It can be used with respiratory protection masks and filters.

Please follow this guide for your safety.



***When using Covidur Face Shield correctly, we recommend that you breathe as above in images.**